



**Vitrotest**<sup>®</sup>

# **Vitrotest TOXO-IgG Avidity**

**ELISA test-kit for the determination of avidity index  
of IgG antibodies to *Toxoplasma gondii***

## **Instruction for use**

1. Intended use
2. Clinical value
3. Principle of the test
4. Materials and equipment
5. Reservations and safety precautions
6. Storage and stability
7. Specimen collection
8. Dilution of samples and reagent preparation
9. Assay procedure
10. Calculation and interpretation of the results
11. Limits of the test

**Reference sources**

**Legend**

IVD

For in vitro diagnostic use

REF

TK064

## «Vitrotest TOXO-IgG Avidity»

### ELISA test-kit for the determination of avidity index of IgG antibodies to *Toxoplasma gondii*

#### 1. Intended use

ELISA test-kit «Vitrotest TOXO-IgG Avidity» is intended for the determination of avidity index of IgG antibodies to *Toxoplasma gondii* in human serum or plasma.

To determine the avidity index in the test «Vitrotest TOXO-IgG Avidity» can only be used sera containing IgG antibodies to *T.gondii* (concentration of specific IgG above 30 IU/ml in the test «Vitrotest TOXO-IgG»).

#### 2. Clinical value

Toxoplasmosis is a parasitic disease, agent of which is an intracellular parasite *Toxoplasma gondii*. *T.gondii* infection rarely leads to clinically significant toxoplasmosis. In case of illness toxoplasmosis occurs latent or chronically. The problem of timely detection of the disease in pregnant women and in patients with suppressed immune systems is of particular importance. The consequences of maternal primary infection during pregnancy in 30-50% of cases lead to abortion or birth of a child with physical and mental disabilities.

At toxoplasmosis specific antibodies of IgM class to *T.gondii* appear on 2nd week of infection, specific IgG antibodies appear later and reach a high concentration at 7-10th week; further high level of IgG keeps or years or reduced to minimum. Increase of more than half titers of IgG class antibodies in serum obtained 2-4 weeks after the initial examination, may indicate either the primary infection or reactivation of chronic toxoplasmosis.

Given the serious consequences of primary toxoplasmosis for pregnancy raises the question of differentiation of primary infection from reinfection or reactivation of toxoplasmosis. Simultaneous detection of specific antibodies of class IgM and IgG does not allow accurate diagnosis of primary toxoplasmosis because specific IgM class antibodies can be detected in human serum for 2 years after infection, as well as reinfection and reactivation disease. In such cases, determination of avidity index of specific IgG antibodies by enzyme immunoassay helps distinguish primary and paste infection. This definition is based on the properties of the immune system change low-avidity synthesis of antibodies during primary infection on the synthesis of high-avidity antibodies at past infection. Under avidity must be understood strength of binding of antigen to the antibody. At early acute primary infection are produced low-avidity IgG antibodies to *T.gondii*. Further, for several months avidity of antibodies increases. High-avidity IgG to *T.gondii* are kept in the body for a long time and detection of them is an evidence of past-infection or reactivation of toxoplasmosis (possibly with simultaneous detection of IgM to *T.gondii*).

#### 3. Principle of the test

Principle of the test of «Vitrotest TOXO-IgG Avidity» kit is based on an indirect enzyme immunoassay technique.

The solid phase is made of strip microplate coated with the purified antigens of *Toxoplasma gondii*. During incubation of samples in wells of ELISA plate specific to *T.gondii* antibodies are bound to the antigen on the solid phase. After aspiration into parallel wells are added dissociation buffer and comparison solution and incubated for a period of time. Denaturing agent in the composition of buffer dissociation destroys the weak links of low-avidity antibodies with antigen. After washing anti-specific conjugate of anti-IgG monoclonal antibodies with horseradish peroxidase added to the wells, it binds to immune complexes in the solid phase. Unbound components are washed out. The immune complex formed in the wells are visualized by adding chromogen solution of 3,3',5,5'- tetramethylbenzidine (TMB). As a result of the reaction solution in wells where immune complexes were formed would be painted in blue. The reaction is stopped by adding stop reagent, blue colored wells become yellow. The result of the analysis is determined by spectrophotometric reading at 450/620 nm. Avidity index is calculated by determining the ratio of optical density in wells with dissociation buffer to an optical density in the wells of comparison solution. The result is expressed as a percentage.

#### 4. Materials and equipment

##### 4.1 Contents of the kit

**ELISA plate** – 12 strips of 8 wells (with the possibility of separation of the wells) with immobilized antigen of *Toxoplasma gondii*.

**High-avidity positive control K+HA** – 1 vial containing 0,5 ml solution of high-avidity specific to *Toxoplasma gondii* human IgG (pink).

**Low-avidity positive control K+LA** – 1 vial containing 0,5 ml solution of low-avidity specific to *Toxoplasma gondii* human IgG (green).

**Negative control** – 1 vial containing 0,5 ml negative human serum (yellow)..

**Sample pre-diluent** – 1 bottle containing 12 ml buffer with detergent and preservatives (brown-green).

**Sample diluent** – 1 bottle containing 12 ml buffer with detergent and preservatives (yellow).

**Conjugate solution** – 1 bottle containing 12 ml buffer solution of monoclonal antibodies to human IgG conjugated with horseradish peroxidase, with stabilizers and preservatives (violet). Ready to use.

**Washing solution Tw20 (20x)** – 1 bottle containing 50 ml 20-fold concentrated phosphate buffer with Tween 20 (colourless).

**Dissociation buffer (DB)** – 1 vial containing 8 ml of denaturing agent. Ready to use (pink).

**Comparison solution (CS)** – 1 vial containing 8 ml of buffer with detergent. Ready to use (blue).

**TMB Solution** – 1 bottle containing 12 ml of TMB solution and hydrogen peroxide, with stabilizers and preservatives (colourless).

**Stop-reagent** – 1 vial containing 12 ml of 0,5M sulphuric acid solution (colourless).

**Plate for pre-dilution of serum** – 12 strips of 8 wells.

**Adhesive film** – 2 sheets of adhesive film for covering the plates during incubation.

**Sera identification plan** – 1 sheet of paper for noting the schemes of samples distribution.

**Instruction for use** – one copy of user manual.

#### 4.2 Additional reagents, materials and equipment

In order to conduct the analysis, the following additional reagents, materials and equipment are required:

- deionized or distilled water;
  - filter paper;
  - graduated cylinders of 10, 200 and 1000 ml capacity;
  - disposable gloves;
  - hydrogen peroxide solution 6%;
  - disposable glassware for preparing the reagents (bottles and trough);
  - timer;
  - mono- and multi-channel automatic adjustable pipettes capable of delivering volumes of 20, 200 and 1000 microliters and tips for them;
  - thermostat for 37 °C;
  - container for solid waste;
  - container for liquid waste;
  - <sup>1</sup>automatic or semi-automatic washer;
  - <sup>2</sup>mono- or multi-channel reader for microplates at 450/620-695 nm.
- <sup>1,2</sup>Please, consult us about the adaptation of your equipment.

## 5. Reservations and safety precautions

### 5.1. Reservations:

- do not use expired reagents;
  - do not use in the analysis and do not mix components of different lots and components of test kits with different nosology;
  - do not use reagents of other manufacturers in composition with the Vitrotest® sets;
  - *Note: possible to use washing solution Tw20 (20X), Sample pre-diluent, TMB solution and Stop-reagent with other series that are different from those attached to the test kit. These reagents are used in other test systems of Ramintek IPC. Please consult us for details.*
  - close reagent vials after use only with their appropriate caps;
  - strictly follow to the washing procedure requirements described in this instruction;
  - control the filling and full aspiration of the solution in the wells;
  - use a new distribution tip for each serum and reagent;
  - protect kit reagents from straight sun rays;
  - TMB solution must be colourless or light blue before its use. If solution is dark blue or yellow, it cannot be used.
- Avoid any contact of the TMB solution with metals or metal's ions. Use glassware thoroughly washed and rinsed with deionized water.
- use only disposable pipette tips for adding TMB-substrate into plate's wells;
  - never use the same glassware for conjugate solution and chromogen.

### 5.2. Safety precautions:

- all reagents included in the kit are intended for "in vitro" diagnostic use;
- wear disposable gloves when perform analysis;
- do not pipette by mouth;
- the controls of «Vitrotest TOXO-IgG Avidity» are negative for anti-HCV, anti-HIV1/2, anti-*T.pallidum* antibodies and HBsAg. Nevertheless, all controls and sera should still be regarded and handled as potentially infectious;
- the liquid waste must be inactivated, for example, with the hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with the sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other disinfectant agents;
- the solid waste must be inactivated with autoclaving at 121°C for 1 hour;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- avoid spilling of TMB-solution and Stop-reagent and any contact of these solutions with the skin or mucosa;
- in case of spilling of solutions, that do not contain acid, e.g. sera, rinse the surface with hydrogen 6% solution, then dry with filter paper.

## 6. Storage and stability

Reagents of the kit are stable up to the expiry date on the label, when stored at 2-8 °C.

Transport the test-kit at 2-8 °C. Disposable transportation at temperature not higher than 20°C during two days is allowed.

## 7. Specimen collection

The serum or plasma samples should be stored at 2-8 °C not more than 3 days after collection. It is possible to store them longer, but frozen (-20 to -70 °C). Before use frozen samples, wait for 30 minutes for the reagents to stabilize at room temperature. Mix thawed samples well to homogeneity. Avoid repeated freezing/thawing. Samples containing aggregates must be clarified by centrifugation. Do not use samples with contaminated, hyperlipemic and hyperhaemolysed sera.

## 8. Dilution of samples and reagent preparation

It is very important to bring all reagents of the «Vitrotest TOXO-IgG Avidity» kit to room temperature (18-25°C) for 30 minutes before use.

### 8.1. ELISA plate preparation

Before opening the ELISA plate, allow it to stabilize at room temperature for 30 minutes to avoid water condensation inside the wells. Open the vacuum bag and remove necessary amount of wells. Immediately after removal of wells, the remaining strips should be **resealed with zip-lock and stored at 2-8 °C**. Microplate in thus packed bag is stable for 3 month.

### 8.2. Washing solution

The vial contains 50 ml of a concentrated 20X buffer with detergent. Dilute the washing solution 1:20 (1+19) with distilled or deionised water, then mix. For example: for 4 ml of concentrate – 76 ml of distilled water is enough for one strip. Crystals in the solution disappear by warming up to 37°C for 15-20 min.

The diluted washing solution can be stored at 2-8°C not more than 7 days.

### 8.3. Pre-dilution of samples and controls

Pre-dilute samples and controls by 10 times with sample pre-diluent. For this purpose, in the required number of wells of plate for pre-dilution of sera (comes in a set) dispense 90 µl of sample pre-diluent and add 10 µl of sera and controls. Gently pipette the mix in wells, avoiding foaming. After addition of serum color of the solution in wells changes from brown-green to blue.

The procedure of dilution of samples and controls should be carried out immediately prior to analysis.

**Attention!** For correct determination of avidity index sera with high level of IgG to *T.gondii* (optical density value in the test «Vitrotest TOXO-IgG» above 3,0 optical units) should be prepared their pre-dilution 1:40. The final dilution of such sera in wells of ELISA plate will be 1: 400.

### 8.4. Dissociation buffer

Dissociation buffer is ready for use. In the case of appearance of crystals in the DB heat the bottle at 37 °C for 15-20 minutes.

## 9. Assay procedure

9.1. Take out from the protective packing the support frame and the necessary number of wells (two wells for each sample and controls). Wells with controls must be included in each analysis.

The figure shows the sequence of dispensation of controls and samples. Thus, the maximum number of samples of serum and/or plasma, which can be studied in the test «Vitrotest TOXO-IgG Avidity» is of 45 analyzes.

	1	2	3	4	5	6	7	8	9	10	11	12
A	K+ HA	K+ HA	6	6	14	14	22	22	30	30	38	38
B	K+ LA	K+ LA	7	7	15	15	23	23	31	31	39	39
C	K-	K-	8	8	16	16	24	24	32	32	40	40
D	1	1	9	9	17	17	25	25	33	33	41	41
E	2	2	10	10	18	18	26	26	34	34	42	42
F	3	3	11	11	19	19	27	27	35	35	43	43
G	4	4	12	12	20	20	28	28	36	36	44	44
H	5	5	13	13	21	21	29	29	37	37	45	45
	CS	DB	CS	DB	CS	DB	CS	DB	CS	DB	CS	DB

9.2. Complete the sera identification plan.

9.3. Prepare washing solution (refer to point 8.2).

9.4. Conduct pre-dilution of sera (refer to point 8.3).

9.5. Dispense 90 µl of sample diluent in each well.

9.6. Dispense into parallel wells of ELISA plate pre-diluted controls and examined samples in the following order: into the wells A1, A2 – 10 µl of pre-diluted 1:10 K+ HA, B1, B2 – 10 µl of pre-diluted 1:10 K+ LA, C1, C2 – 10 µl of pre-diluted 1:10 negative control, D1, D2 – 10 µl of pre-diluted 1:10 sample №1 etc., according to the scheme. Thus, the final dilution of serum in wells of ELISA plate should be 1: 100. Gently pipette the mix in wells, avoiding foaming – color of the solution in wells changes from yellow to green.

9.7. Cover strips with an adhesive film and incubate for 30 min at 37°C.

9.8. After completing the incubation remove the adhesive film carefully and aspirate the content of wells using the automatic washer or 8-channel pipette.

9.9. Disperse 100 µl of comparison solution to the wells of odd stripes and 100 µl of dissociation buffer to the wells of even stripes. Cover strips with an adhesive film, incubate for 15 min at room temperature (18-25 °C).

9.10. After completing the incubation remove the adhesive film carefully and wash the wells five times using the automatic washer or 8-channel pipette in the following order:

– aspirate the content of wells strips into a liquid waste container;

– fill the strip wells with a minimum of 300 microliters of washing solution for each well (respect the soak time of a minimum of 30 seconds);

– aspirate the solution of all wells, the residual volume of solution after aspiration must be lower than 5 microliters;

– repeat the washing step four more times;

– after the last aspiration blot the microplate by turning it upside down on absorbent paper.

9.11. Dispense 100 µl of conjugate solution per well (into all wells). Cover strips with an adhesive film, incubate for 30 min at 37°C.

9.12. After completing the incubation remove the adhesive film carefully and wash the wells five times as described above (refer to point 9.10).

9.13. TMB-solution is ready to use TMB-substrate solution with hydrogen peroxide. TMB-solution should be colorless, protect TMB-substrate solution from straight sun rays. To add TMB-solution only new tips must be used: carefully select a TMB-solution from the vial and without touching the bottom and walls of the hole plate, add 100 µl TMB solution per well.

9.14. Incubate the strips for 30 minutes at room temperature of 18-25°C in the dark. Do not use adhesive film in this incubation.

9.15. Add 100 µl of stop-reagent in each well. Respect the same distribution sequence than for the TMB-substrate solution.

9.16. Read the optical density of every strip well in dual wavelength reading at 450/620 nm, within the 5 minutes of stopping the reaction. Pay attention to the cleanness of the wells bottom outside.

*Measurement in the single-wave procedure at 450 nm is possible. Reserve blank well to adjust spectrophotometer in such analysis. Only TMB-substrate and stop-reagent must be added in blank well.*

## 10. Calculation and interpretation of the results

### 10.1. Test validity conditions:

In order for an assay to be considered valid, the following criteria must be met:

– optical density (OD) of HA positive control in the well treated by CS is not lower than 1,5 optical units (OU),

– OD of LA positive control in the well treated by CS is not lower than 0,8 OU,

– OD of negative control should be lower than 0,15 OU,

– avidity index of K+ HA > 40%

– avidity index of K+ LA < 30%

### 10.2. Calculation of the results.

Avidity index is calculated as:

$$AI = \frac{OD\ DB}{OD\ CS} \times 100\%$$

where AI – avidity index, OD DB – optical density of the sample in the well with dissociation buffer, OD CS – optical density of the sample in the well with comparison solution.

### 10.3. Interpretation of the results.

For sera with OD value below 0,300 OU in the well with comparison solution avidity index is not calculated.

Results of determination of avidity index of specific IgG to *T.gondii* are interpreted as follows:

<i>AI</i>	<i>Result</i>
< 30 %	Sample contains low-avidity IgG
30-40%	Indeterminate result
> 40 %	Sample contains high-avidity IgG



### **11. Limits of the test**

Detection of low-avidity IgG in the test «Vitrotest TOXO-IgG Avidity» is an evidence of primary infection *T.gondii*.

Detection of high-avidity IgG in the test «Vitrotest TOXO-IgG Avidity» is an evidence of past-infection or reinfection *T.gondii*.

Upon indeterminate results should re-examine the sample and in case of indeterminate result again explore new specimen obtained from this patient as soon as possible.

For diagnosis should be considered as the results of laboratory tests and clinical manifestations of the disease.

### **Reference sources**

1. Дубинська Г.М., Ізюмська О.М., Козюк П.М. та ін. //Токсоплазмоз. Клініка, діагностика та лікування. // Вести медицины и фармации. – 2003. – №1. – С. 23-25.

2. Desmonts G., Daffos F., Forestier F. et al. Prenatal Diagnosis of Congenital Toxoplasmosis // Lancet. – 1985. – Vol. 8427, N1. – P. 500-504.












3. Genum P.A., Stray-Pedersen B., Gundersen A. Improved Diagnosis of Primary Toxoplasma Gondii Infection in Early Pregnancy by Determination Antitoxoplasma Immunoglobulin G Avidity // Journal of Clinical Microbiology. – 1997. – N8 – P. 1972-1977.

4. S.Singh. Mother-to-child transmission and diagnosis of Toxoplasma Gondii infection during pregnancy // Indian Journal of Medical Microbiology. – 2003. – Vol. 21, N2. – P.69-76.

5. Tenter A.M., Heckeroth A.R., Weiss L.M. Toxoplasma gondii: from animals to humans // International Journal for Parasitology. – 2000. – №30 – P. 1217-1258.

## Legend

Interpretation of notation conventions:

	For in vitro diagnostic use
	Batch code
	Catalogue number
	Production date
	Expiry date
	Storage temperature limitation
	Meant for <n> tests
	Protect from direct solar radiation
	Attention! Consult instruction for use
	Manufacturer and its address
	Consult instructions for use

*For questions and suggestions regarding the kit, contact the manufacturer:*



Ramintek Innovation-Production Company  
03039 Ukraine, Kiev, 7 October 40th Anniversary Av., of. 227 (registered address)  
07300 Vishgorod, Kiev region, 19 Sholudenko Str. (factual address)  
Tel. +380 44 222-76-72

## Scheme of the assay «Vitrotest TOXO-IgG Avidity»

---

Keep reagents at room temperature 18-25°C during 30 minutes

---

Prepare washing solution, dilute 20x concentrate washing solution *Tw20* with distilled water 1:20 (1+19).  
For example, 4 ml of solution + 76 ml of water

---

Complete the sera identification plan

---

Prepare pre-dilution (1:10) of sera: in wells of plates for pre-dilution add 90 µl of sample pre-diluent (brown-green) and 10 µl of samples or controls.

*After dispensing of serum the color in well switches from brown-green to blue*

---

Dispense 90 µl of sample diluent (yellow) and add 10 µl of pre-diluted controls and patient samples into the parallel wells:

A1, A2 – positive control HA,    B1, B2 – positive control LA,  
C1, C2 – negative control,    D1, D2 etc. – patient samples

---

Cover wells with an adhesive film, incubate for 30 min at 37°C

---

Aspirate content of wells

---

Into wells of odd strips dispense 100 µl of comparison solution (blue) and into wells of even strips – 100 µl of dissociation buffer (pink)

---

Cover wells with an adhesive film, incubate for 15 min at 18-25°C

---

Wash wells five times

---

Dispense 100 µl of conjugate solution (violet) into the wells

---

Cover wells with adhesive film, incubate for 30 min at 37°C

---

Wash wells five times

---

Dispense 100 µl of TMB substrate solution into the wells

---

Incubate the plate for 30 min in the dark at room temperature (18-25°C)

---

Add 100 µl of stop-reagent in each well

---

Read optical density at 450/620 nm

---

Calculate avidity index of IgG to *Toxoplasma gondii* for samples

---

Interpret the results according to the table

<i>AI</i>	<i>Result</i>
< 30 %	Sample contains low-avidity IgG
30-40%	Indeterminate result
> 40 %	Sample contains high-avidity IgG