

The logo for Vitrotest, featuring a stylized 'Y' shape on the left and the word 'Vitrotest' in a bold, sans-serif font to its right, with a registered trademark symbol (®) to the upper right of the 't'.

**Vitrotest**<sup>®</sup>

# **Vitrotest Anti-Trichinella**

**ELISA test-kit for the detection of antibodies specific to *Trichinella spiralis***

## **Instruction for use**

1. Intended use
2. Clinical value
3. Principle of the test
4. Materials and equipment
5. Reservations and safety precautions
6. Storage and stability
7. Specimen collection
8. Reagent preparation
9. Assay procedure
10. Calculation and interpretation of the results
11. Performance characteristics
12. Limits of the test

**Reference sources**

**Legend**

IVD

For in vitro diagnostic use

REF

TK067

## «Vitrotest Anti-Trichinella»

### ELISA test-kit for the detection of antibodies specific to *Trichinella spiralis*

#### 1. Intended use

ELISA test-kit «Vitrotest Anti-Trichinella» is intended for the detection of antibodies IgG and IgA classes to *Trichinella spiralis* in human serum or plasma.

The test kit may be used both for the ELISA using automatic pipettes and standard equipment and for setting with the open-system immunoenzymatic automated analyzer.

#### 2. Clinical value

Trichinosis is a helminthiasis, which is caused by nematode *Trichinella spiralis*. *Trichinella* are small, almost thread-like worms that are covered with cross-striped cuticle. Length of sexually mature male is 1,2-2 mm, thickness - 0,04-0,05 mm. Length of mature females before fertilize is 1,5-1,8 mm, after fertilization it increases to 4,4 mm.

Infection occurs by ingestion of meat containing live encapsulated trichinella larvae. In the process of digestion by the activity of gastric juice larvae released from capsules to penetrate the submucosal layer of the small intestine, fixed on the mucosa and begin to proliferate. After fertilization, the males die, and females 2 days after infestation begin to bear larvae, which enter the blood and lymphatic vessels through tissue mucosa and spread throughout the body and deposit in striated muscle. After 18-20 days after infection larvae in muscles lengthened to 0,8 mm, reach invasive stage and begin to convolve in the form of a spiral. As an answer of surrounding tissues in 35-40 days around the larvae is formed connective-tissue capsule, which later is impregnated with calcium salts that leads to calcification. The larvae remain viable for many years.

The incubation period for human trichinosis lasts 10-25 days. *Trichinella* initially is localized in the intestine. Within 1-2 days after infection there are symptoms such as nausea, heartburn, indigestion, diarrhea, severity of which depends on the degree of damage. Further symptoms and the degree of their manifestation depend on the location of parasite in different parts of the body. Trichinosis is characterized by fever, myalgia, facial swelling, skin rash, blood eosinophilia, and in severe leakage - damage internal organs and the central nervous system.

Trichinosis is not always immediately recognizable, especially in areas where is rare. Trichinosis diagnosis is based on clinical, epidemiological history, serological tests, such as complement fixation, indirect hemagglutination reaction and ELISA. The latter method is recommended by OIE for the serological diagnosis of trichinosis.

The most specific method to confirm the infestation is to detect IgG antibodies to trichinella antigens in the blood, which can be determined from 2-3 to 4-6 weeks after eating contaminated meat. Specific IgE class antibodies are present in the blood at the acute stage of the disease, however, are rare, due to the short period of their circulation in the bloodstream. With early clinical manifestations of invasion specific antibodies might not be detected. Therefore, the suspected trichinosis and a negative result the research is carried out again in a week or two. Seroconversion usually occurs during the second-fifth week of infections, depending on the infectious dose received by the patient. Assessment of the dynamics of antibody is very informative criterion of effectiveness of the therapy - with ineffective or deferred therapy specific antibodies are detected for years (up to 20 years); if effective treatment in the first two weeks after infection, antibodies disappear within a year.

#### 3. Principle of the test

Principle of the test of «Vitrotest Anti-Trichinella» kit is based on an indirect enzyme immunoassay technique.

The solid phase is made of strip microplate coated with antigens of *Trichinella spiralis* larvae. During incubation of samples in wells of ELISA plate specific to *Trichinella spiralis* antibodies are bound to the antigen on the solid phase. After washing out unbound components the mixture of anti-specific conjugates of anti-IgG and anti-IgA monoclonal antibodies with horseradish peroxidase added to the wells, it binds to immune complexes in the solid phase. Unbound components are washed out. The immune complex formed in the wells are visualized by adding chromogen solution of 3,3',5,5'- tetramethylbenzidine (TMB). As a result of the reaction solution in wells where immune complexes were formed would be painted in blue. The reaction is stopped by adding stop reagent, blue colored wells become yellow. The result of the analysis is determined by spectrophotometric reading at 450/620 nm.

#### 4. Materials and equipment

##### 4.1 Contents of the kit

**ELISA plate** – 12 strips of 8 wells (with the possibility of separation of the wells) with immobilized antigens of *Trichinella spiralis* larvae.

**Positive control** – 1 vial containing 0,3 ml solution of human antibodies specific to *Trichinella spiralis* (pink).

**Negative control** – 1 vial containing 0,5 ml negative human serum (yellow).

**Sample diluent** – 1 bottle containing 12 ml buffer with skim milk extract, detergent and preservatives (brown-green).

**Conjugate solution** – 1 bottle containing 12 ml buffer solution of monoclonal antibodies to human IgG conjugated with horseradish peroxidase, with stabilizers and preservatives (green). Ready to use.

**TMB Solution** – 1 bottle containing 12 ml of TMB solution and hydrogen peroxide, with stabilizers and preservatives (colourless).

**Washing solution Tw20 (20x)** – 1 bottle containing 50 ml 20-fold concentrated phosphate buffer with Tween-20 (colourless).

**Stop-reagent** – 1 vial containing 12 ml of 0,5 M sulphuric acid solution (colourless).

**Adhesive film** – 2 sheets of adhesive film for covering the plates during incubation.

**Sera identification plan** – 1 sheet of paper for noting the schemes of samples distribution.

**Instruction for use** – one copy of user manual.

#### 4.2 Additional reagents, materials and equipment

In order to conduct the analysis, the following additional reagents, materials and equipment are required:

- deionized or distilled water;
  - filter paper;
  - graduated cylinders of 10, 200 and 1000 ml capacity;
  - disposable gloves;
  - hydrogen peroxide solution 6%;
  - disposable glassware for preparing the reagents (bottles and trough);
  - timer;
  - mono- and multi-channel automatic adjustable pipettes capable of delivering volumes of 20, 200 and 1000 microliters and tips for them;
  - thermostat for 37 °C;
  - container for solid waste;
  - container for liquid waste;
  - <sup>1</sup>automatic or semi-automatic washer;
  - <sup>2</sup>mono- or multi-channel reader for microplates at 450/620-695 nm.
- <sup>1,2</sup>Please, consult us about the adaptation of your equipment.

## 5. Reservations and safety precautions

### 5.1. Reservations:

- do not use expired reagents;
  - do not use in the analysis and do not mix components of different lots and components of test kits with different nosology;
  - do not use reagents of other manufacturers in composition with the Vitrotest® sets;
  - *Note: possible to use washing solution Tw20 (20X), TMB solution and Stop-reagent with other series that are different from those attached to the test kit. These reagents are used in other test systems of Ramintek IPC. Please consult us for details.*
  - close reagent vials after use only with their appropriate caps;
  - strictly follow to the washing procedure requirements described in this instruction;
  - control the filling and full aspiration of the solution in the wells;
  - use a new distribution tip for each serum and reagent;
  - protect kit reagents from straight sun rays;
  - TMB solution must be colourless or light blue before its use. If solution is dark blue or yellow, it cannot be used.
- Avoid any contact of the TMB solution with metals or metal's ions. Use glassware thoroughly washed and rinsed with deionized water.
- use only disposable pipette tips for adding TMB-substrate into plate's wells;
  - never use the same glassware for conjugate solution and chromogen.

### 5.2. Safety precautions:

- all reagents included in the kit are intended for "in vitro" diagnostic use;
- wear disposable gloves when perform analysis;
- do not pipette by mouth;
- the controls of «Vitrotest Anti-Trichinella» are negative for anti-HCV, anti-HIV1/2, anti-*T.pallidum* antibodies and HBsAg. Nevertheless, all controls and sera should still be regarded and handled as potentially infectious;
- the liquid waste must be inactivated, for example, with the hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with the sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other disinfectant agents;
- the solid waste must be inactivated with autoclaving at 121°C for 1 hour;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- avoid spilling of TMB-solution and Stop-reagent and any contact of these solutions with the skin or mucosa;
- in case of spilling of solutions, that do not contain acid, e.g. sera, rinse the surface with hydrogen 6% solution, then dry with filter paper.

## 6. Storage and stability

Reagents of the kit are stable up to the expiry date on the label, when stored at 2-8 °C.

Transport the test-kit at 2-8 °C. Disposable transportation at temperature not higher than 20°C during two days is allowed.

## 7. Specimen collection

The serum or plasma samples should be stored at 2-8 °C not more than 3 days after collection. It is possible to store them longer, but frozen (-20 to -70 °C). Before use frozen samples, wait for 30 minutes for the reagents to stabilize at room temperature. Mix thawed samples well to homogeneity. Avoid repeated freezing/thawing. Samples containing aggregates must be clarified by centrifugation. Do not use samples with contaminated, hyperlipemic and hyperhaemolysed sera.

## 8. Reagent preparation

*It is very important to bring all reagents of the «Vitrotest Anti-Trichinella» kit to room temperature (18-25°C) for 30 minutes before use.*

### 8.1. ELISA plate preparation

Before opening the ELISA plate, allow it to stabilize at room temperature for 30 minutes to avoid water condensation inside the wells. Open the vacuum bag and remove necessary amount of wells. Immediately after removal of wells, the remaining strips should be **resealed with zip-lock and stored at 2-8 °C**. Microplate in thus packed bag is stable for 3 month.

### 8.2. Washing solution

The vial contains 50 ml of a concentrated 20X buffer with detergent. Dilute the washing solution 1:20 (1+19) with distilled or deionised water, then mix. For example: for 4 ml of concentrate – 76 ml of distilled water is enough for one strip. Crystals in the solution disappear by warming up to 37°C for 15-20 min.

The diluted washing solution can be stored at 2-8°C not more than 7 days.

## 9. Assay procedure

9.1. Take out from the protective packing the support frame and the necessary number of wells (the number of investigated samples and four wells for controls). Wells with the controls must be included in each analysis.

9.2. Complete the sera identification plan.

9.3. Prepare washing solution (refer to point 8.2).

9.4. Dispense 90 µl of sample diluent in each well.

9.5. Dispense 10 µl of controls and unknown samples into the wells in the following order: A1 – positive control, B1, C1 and D1 – negative control, E1 and rest wells – unknown samples. Gently pipette the mix in wells, avoiding foaming. After addition of serum color of the solution in wells changes from brown-green to blue.

9.6. Cover strips with an adhesive film and incubate for 30 min at 37 °C.

9.7. After completing the incubation remove the adhesive film carefully and wash the wells five times using the automatic washer or 8-channel pipette in the following order:

– aspirate the content of wells strips into a liquid waste container;

– fill the strip wells with a minimum of 300 microliters of washing solution for each well (respect the soak time of a minimum of 30 seconds);

– aspirate the solution of all wells, the residual volume of solution after aspiration must be lower than 5 microliters;

– repeat the washing step four more times;

– after the last aspiration blot the microplate by turning it upside down on absorbent paper.

9.8. Dispense 100 µl of conjugate solution per well. Cover strips with an adhesive film, incubate for 30 min at 37 °C.

9.9. After completing the incubation remove the adhesive film carefully and wash the wells five times as described above (refer to point 9.7).

9.10. TMB-solution is ready to use TMB-substrate solution with hydrogen peroxide. TMB-solution should be colorless, protect TMB-substrate solution from straight sun rays. To add TMB-solution only new tips must be used: carefully select a TMB-solution from the vial and without touching the bottom and walls of the hole plate, add 100 µl TMB solution per well.

9.11. Incubate the strips for 30 minutes at room temperature of 18-25°C in the dark. Do not use adhesive film in this incubation.

9.12. Add 100 µl of stop-reagent in each well. Respect the same distribution sequence than for the TMB-substrate solution.

9.13. Read the optical density of every strip well in dual wavelength reading at 450/620 nm, within the 5 minutes of stopping the reaction. Pay attention to the cleanness of the wells bottom outside.

*Measurement in the single-wave procedure at 450 nm is possible. Reserve blank well to adjust spectrophotometer in such analysis. Only TMB-substrate and stop-reagent must be added in blank well.*

## 10. Calculation and interpretation of the results

### 10.1. Test validity conditions:

Calculate the mean optical density (OD) of negative control

$$OD_{NC_{mean}} = (OD_{NC1} + OD_{NC2} + OD_{NC3})/3.$$

In order for an assay to be considered valid, the following criteria must be met:

- OD of the positive control is not lower than 0,8 optical unit (OU),
- OD of negative controls should be lower or equal to 0,15 OU,
- OD of every negative control differs no more than twice from the mean value of negative control that is

$$OD NC_{mean} \times 0,5 \leq OD NC_n \leq OD NC_{mean} \times 2,0.$$

If one of the negative controls does not respect this norm, disregard and recalculate the mean using remaining values.

#### 10.2. Calculation of the results.

Calculate cut-off by adding value 0,30 to the mean NC, that is

$$Cut\ off = OD\ NC_{mean} + 0,30.$$

Calculate the index of positivity (IP) for each sample

$$IP = \frac{OD_{sample}}{Cut\ off}$$

#### 10.3. Interpretation of the results

The samples with IP above 1,1 are considered **positive (IP > 1,1)**.

The samples with IP below 0,9 are considered **negative (IP < 0,9)**.

The samples with IP **within 0,9-1,1** are considered **indeterminate (0,9 ≤ IP ≤ 1,1)**. It is recommended to retest such samples in duplicate. If the results are again within indeterminate, it is necessary to test sera obtained after 2-4 weeks. If you get undefined results assume that the serum does not contain specific antibodies to *Trichinella spiralis*.

## 11. Performance characteristics

### 11.1. Specificity and sensitivity

Sensitivity of the test «Vitrotest Anti-Trichinella» was evaluated by using panel of characterized sera, consisting of 71 samples of human blood sera containing antibodies to *Trichinella spiralis* (sera were pre-tested in a different commercial test kit for the detection of IgG antibodies to *Trichinella spiralis*). In the test «Vitrotest Anti-Trichinella» 67 sera were identified as positive. In the study of 152 negative for antibodies to *Trichinella spiralis* sera specificity rate was over 98%.

### 11.2. Accuracy

#### Intra assay reproducibility

Coefficient of variation (CV) for two sera with different levels of specific antibodies was calculated in 32 repetitions on two series of test kits.

Serum №	IP	CV <sub>1</sub> , %	CV <sub>2</sub> , %
009	2,1	4,3	3,2
034	6,5	5,4	4,5

#### Inter assay reproducibility

Coefficient of variation (CV) for two sera with different levels of specific antibodies was calculated for three days in three ELISA performances, in eight repetitions for each analysis.

Serum №	IP	CV, %
009	2,0	3,7
034	6,6	5,1

## 12. Limits of the test

A positive result in the test kit «Vitrotest Anti-Trichinella» is an indication of presence in patient of antibodies of class IgG, specific to *Ascaris lumbricoides*. The presence of antibodies of this class in infants is not evidence of invasion of *Trichinella spiralis*.

Undefined results may indicate *Trichinella spiralis* infestation in history.

A negative result in the «Vitrotest Anti-Trichinella» kit indicates the absence of antibodies to *Trichinella spiralis* in studied human serum, or the concentration of specific antibodies is below the sensitivity analysis. Specific antibodies cannot be detected in early clinical manifestations of infestation. In this case it is recommended in a week or two get and check serum samples from the patients with clinical signs of trichinosis.












The final diagnosis cannot be established only on the basis of serological test. At diagnosis should take into account the complex of laboratory and instrumental studies, and clinical manifestations of the disease. It is impossible to completely eliminate cross-reactions with antibodies to antigens of other helminths. To exclude false-positive results is recommended verification study of positive samples by immunoblotting.

## Reference

1. Гельминтозы человека // Медицинская микробиология / Гл. ред. В.И. Покровский, О.К.Позеев – М.: ГЭОТАР Медицина, 1998. – С. 942-986.
2. Руководство по инфекционным болезням с атласом инфекционной патологии. Под редакцией Ю.В. Лобзина, С.С. Козлова, А.Н. Ускова. – СПб, 2000.
5. Dupouy-Camet J. Trichinellosis: a worldwide zoonosis. // *Vet Parasit*, - 2000, - V.93 - p.191-200.
6. Dupouy-Camet, J. and Bruschi, F. Management and diagnosis of human trichinellosis. In, (Dupouy-Camet, J and Murrell, K.D. eds.), *FAO/WHO/OIE Guidelines for the Surveillance, Management, Prevention and Control of Trichinellosis*, Rome, - 2007. - p. 37-68.
7. Gamble HR, Pozio E, Bruschi F, Nockler K, Kapel CM, Gajadhar AA. International Commission on Trichinellosis: recommendations on the use of serological tests for the detection of *Trichinella* infection in animals and man. // *Parasite*. – 2004. - Mar; 11(1) – p.3-13.
8. Gomes-Morales M.A., Ludovisi A., Amati A., et. al. Validation of an enzyme-linked Immunosorbent assay for diagnosis of human trichinellosis. // *Clinical and Vaccine Immunology*, - 2008. - Vol.15 No.11 - p.1723-1729.
9. Murrell K.D., & Pozio E. Trichinellosis: the zoonosis that won't go quietly. // *International Journal of Parasitology*, - 2000, - 30 (12-13), - p.1339-49
10. Kotula A.W., Sharar A., Paroczay E., Gamble H.R., Murrell K.D. & Douglass L. Infectivity of *Trichinella* from frozen pork. // *Journal of Food Protection* - 1990. – Vol.53 – p.571-573.

## Legend

Interpretation of notation conventions:

	For in vitro diagnostic use
	Batch code
	Catalogue number
	Production date
	Expiry date
	Storage temperature limitation
	Meant for <n> tests
	Protect from direct solar radiation
	Attention! Consult instruction for use
	Manufacturer and its address
	Consult instructions for use

*For questions and suggestions regarding the kit, contact the manufacturer:*



Ramintek Innovation-Production Company  
03039 Ukraine, Kiev, 7 October 40th Anniversary Av., of. 227 (registered address)  
07300 Vishgorod, Kiev region, 19 Sholudenko Str. (factual address)  
Tel. +380 44 222-76-72



## Scheme of the assay «Vitrotest Anti-Trichinella»

---

Keep reagents at room temperature (18-25°C) during 30 minutes

---

Prepare washing solution, dilute 20x concentrate washing solution *Tw20* with distilled water 1:20 (1+19).  
For example, 4 ml of solution + 76 ml of water

---

Complete the sera identification plan

---

Dispense 90 µl of sample diluents into the wells

---

Dispense 10 µl of controls and patient samples into the wells:

A1 – positive control,

B1, C1, D1 – negative control,

E1 and other wells – patient samples

*After dispensing of serum the color in well switches from brown to blue*

---

Cover wells with adhesive film, incubate for 30 min at 37°C

---

Wash wells five times

---

Dispense 100 µl of conjugate solution (green) into the wells

---

Cover wells with adhesive film, incubate for 30 min at 37°C

---

Wash wells five times

---

Dispense 100 µl of TMB substrate solution into the wells

---

Incubate the plate for 30 min in the dark at room temperature (18-25°C)

---

Add 100 µl of stopping solution in each well

---

Read optical density at 450/620 nm

---

Calculate the cut-off of the assay «Vitrotest Anti-Trichinella»:

$$\text{Cut-off} = OD \text{ NC mean} + 0,3$$

---

Calculate the index of positivity (IP) for patient samples:  $IP = \frac{OD \text{ of patient sample}}{\text{cut off}}$

---

Interpret the results according to the table:

IP value	Result
$IP_{\text{sample}} > 1,1$	positive
$0,9 \leq IP_{\text{sample}} \leq 1,1$	indeterminate
$IP_{\text{sample}} < 0,9$	negative